

Stainless Steel Auxiliary Products

SS Closed and Open Spools



An economical offering, Stainless Steel springs provide high initial forces; however, unlike the forces of our premium Nickel Titanium springs, stainless steel forces of our premium Nickel Titanium springs, stainless steel forces decrease rapidly with tooth movement. Stainless Steel springs are also subject to taking a permanent set when fully compressed or expanded.

40 inch spool

	CLOSED	OPEN
.010" x .028"		ST-SPC 10x28
.010" x .030"	ST-SPE10x30	ST-SPC10x30
.010" x .036"		ST-SPC10x36

- Extension or closed springs are effective in maintaining space openings.
- Compression or open springs are most commonly used to create space.

SS Kobayashi Ties



SS Ligature Ties



Our Preformed Ligature and Kobayashi products are made of annealed 304 Stainless Steel to ensure consistent shape and predictable performance. Used in place of elastic ligatures, you will appreciate the ease with which they can be placed over the bracket and tied.

Stainless Steel Kobayashi Ties

	Long	Short
	100 pieces	500 pieces
.010"	ST-KTL10	
.012"	ST-KTL12	ST-KTS12
.014"	ST-KTL14	ST-KTS14

Stainless Steel Ligature Ties

		Long	Short
1000 pieces	.008"	ST-LWL8	
1000 pieces	.009"	ST-LWL9	
500 pieces	.010"	ST-LWL10	ST-LWS10
1000 pieces	.012"	ST-LWL12	

Wire Accessories

SS Crimpable Ball Hooks



Crimpable Stop	Item #
RIGHT- PK 10	SSBHR
LEFT- PK 10	SSBHL
STRAIGHT- PK 10	SSBHS

Stainless Steel Stop Lock



Description	Item #
STOP LOCK - PK 5	SL
WRENCH	SLW

SS Crimpable C Stops



Description	Item #
C Stop- PK 10	C-STOP

SS Crimpable Split Stops



Description	Item #
SPLIT STOP- PK 10	S