

ARCHWIRE AVAILABILITY

TRIDENT™ CNA Beta III/TM Nickel Free

ARCH FORM	Round	Rectangular	pc/Pack
EURO			10
NATURAL			10
STANDARD			10
UNIVERSAL			10
STRAIGHT 14 inch			10
Mushr. 2-Loops			10
CTA Intrusion Arch			10

Mushroom Looped available for:
 Upper arch form: 32-56mm in 2mm increments
 Lower arch form: 24-30mm in 2mm increments
 Anterior Loop Spacing is measured from center of loops

ARCHWIRE SIZES

UNIVERSAL available only in .016x.022 - .016x.025 - .017x.025 - .018x.025 - .019x.025
 STRAIGHT length wires are available in all of the above wire sizes as well as .027 - .032 - .036
 CTA INTRUSION ARCH available only in .016x.022 and .017x.025



Round	.016	.018			
Square	.016 x .016	.0175 x .0175	.018 x .018		
Square	.016 x .022	.017 x .025	.018 x .025	.019 x .025	.021 x .025

Single Packaged Wires available for: EURO, NATURAL, UNIVERSAL

Single packaged archwires on request available, except looped archwires.

10 individually packaged wires per pack

*** Minimum order of 10 packs per item number required ***



ARCHWIRE ORDER INFORMATION

In order to facilitate your ordering process, you can customize your wire and create your item number by combining the following digits:

CNA + Arch Form + Extra + Wire Diameter + Jaw Code

40°	OBLIGATORY		OPTIONAL		OBLIGATORY		OBLIGATORY	
35°	ARCH FORM		Extra		Wire Diameter		Jaw Code	
CNA	EURO	E	No Extra		.012"	12	Upper	U
	NATURAL	N	24mm	24	.014"	14	Lower	L
	STANDARD	S	26mm	26	.016"	16	Upper & Lower	
	UNIVERSAL	U	28mm	28		
	STRAIGHT	ST019" x .025"	19x25		
	MUSHR. 2 Loop	M	56mm	56				
	CTA INTR. ARCH	CTA						

Please refer to the corresponding Arch-wire Availability to see which arch diameters are available for your selection!

CNA + N + 14 + U = CNAN14U

EXAMPLE 1 You would like to order a **CNA Beta III Nickel Free** archwire in **NATURAL** Form with the diameter **.014** for the **Upper Arch**.

CNA + M + 34 + 19x25 = CNAM3419x25

EXAMPLE 2 You would like to order a **CNA Beta III Nickel Free** archwire with **Mushroom 2 Loop** with **34mm Anterior Loop Spacing** in wire diameter **.019 x .025**.

For Single Packaged Arch Wires please use the same Order Logic as above stated. Please put a "S-" in front of the item number to indicate that you would like to have your wire selection individually packaged.