

LINGUAL ATTACHMENTS sizes and order Info



Lingual Sheaths

Our Lingual Sheaths are designed to accept palatal bars, quad helix appliances, and other maxillary and mandibular removable format appliances incorporating .036" x .072" distal engagement ends. Locking indent and optional soldered post ensure secure retention.

Lingual Weldable Sheaths

		Size	Item #			Item #
	Lingual sheath with no window, with indent	2 x 0.9 mm / 2 x .035"	LGS-UI		Lingual Sheath with window indent, and mesial ball hook	UR/LL Ea. UL/LR Ea. LGS-WIM-1 LGS-WIM-2
	Lingual Sheath with window indent	2 x 0.9 mm / 2 x .035"	LGS-UIW		Lingual Sheath with window indent, and distal ball hook	UR/LL Ea. UL/LR Ea. LGS-WIO-1 LGS-WIO-2

Lingual Button

		Straight Item #	Curved Item #
	Weldable	AD-LGB-G	AD-LGB-R
	Bondable	AD-LGB-GB	AD-LGB-BB

Lingual Weldable Double Cleat & Eyelets

		Item #
	Lingual double cleat weldable large 8.1mm	LGC-L
	Lingual double cleat weldable small 5.3mm	LDH-S
	Weldable Eyelets .036 ID Small Ea.	M005-1LGE-S
	Weldable Eyelets .036 ID Large Ea.	M005-1LGE-L