

## ARCHWIRE AVAILABILITY

### THERMADENT™ HEAT ACTIVATED Nickel Titanium Archwires 35° and 40°

THERMADENT™ is available in two activation ranges:

#### 35°

THERMADENT™ Medium Force is our specially developed heat-activated Nickel Titanium wire with consistent  $A_1$  performance and moderate forces. This is an easy wire to work with being activated at body temperature.

#### 40°

THERMADENT™ Light Force is a heat-activated wire specifically designed for consistent  $A_1$  performance and low tooth moving forces. THERMADENT™ wire provides outstanding resiliency and exhibits true thermal performance at temperatures higher than body temperature.

ARCH FORM	Round	Rectangular	Round-Dimpled	Rectangular-Dimpled	Pre-Stopped Round	Pre-Stopped Rectangular	pc/Pack
EURO	35°	40° + 35°	40° + 35°	40° + 35°	40° + 35°	40° + 35°	10
NATURAL	35°	40° + 35°	40° + 35°	40° + 35°	40° + 35°	40° + 35°	10
STANDARD	35°	35°					10
UNIVERSAL	35°						10
BIO-OVOID	35°	35°					10
STRAIGHT	35°						10
RC 1, 3, 5	40°	40°	40°	40°			10

\*\*\* PRE-Stopped wires are with crimpable tube stops. \*\*\*  
**.020 and .021 x .025 NOT available**

### ARCHWIRE SIZES

**.012 and .020 UNIVERSAL not available**

Round	.012	.014	.016	.018	.020
Square	.016 x .016	.0175 x .0175	.018 x .018	.020 x .020	
Rectangular	.014 x .025	.016 x .022	.016 x .025	.017 x .025	.018 x .025
	.019 x .025	.021 x .025			

**Single Packaged Wires available for:** EURO, NATURAL, BIO-OVOID, UNIVERSAL, RC 1 + 3 + 5

Single packaged archwires on request available.

10 individually packaged wires per pack.

\*\*\* Minimum order of 10 packs per item number required \*\*\*

## ARCHWIRE ORDER INFORMATION

In order to facilitate your ordering process, you can customize your wire and create your item number by combining the following digits - If you would like to order Heat Activated NiTi in 40°, please use NNTM for ordering, if you would like to receive Heat Activated NiTi in 35°, please use NNT.

**NNT or NNTM + Arch Form + Extra + Wire Diameter + Jaw Code**

40° 35°	OBLIGATORY ARCH FORM	OPTIONAL Extra	OBLIGATORY Wire Diameter	OBLIGATORY Jaw Code
NNTM	EURO	E	Tube Stopped	T
	NATURAL	N	Dimpled	D
	STANDARD	S	No Extra	
	UNIVERSAL	U		
	BIO-OVOID	O		
or	STRAIGHT	ST		
NNT	REVERSE CURVE 1	RC1		
	REVERSE CURVE 3	RC3		
	REVERSE CURVE 5	RC5		
			.012"	12
			.014"	14
			.016"	16
			...	...
			.016" x .016"	16x16

Please refer to the table for wire Extras available for your selected arch form!

Please refer to the corresponding Arch-wire Availability to see which arch diameters are available for your selection!

**NNT + E + T + 14 + U = NNTET14U**

**EXAMPLE 1** You would like to order a **NiTi Heat Activated 35°** archwire in **EURO** Form with **Tube Stops** on the archwire with the diameter **.014** for the **Upper Arch**.

$$\text{NNTM} + \text{U} + 16\text{x}16 = \text{NNTMU16x16}$$

**EXAMPLE 2** You would like to order a **NiTi Heat Activated 40°** archwire in **UNIVERSAL** Form with the diameter **.016x.016** for the **Upper and Lower Arch**.

For Single Packaged Arch Wires please use the same Order Logic as above stated. Please put a "S-" in front of the item number to indicate that you would like to have your wire selection individually packaged.