


## ARCHWIRE AVAILABILITY

### TWISTADENT™ Stainless Steel Twisted/Braided/Stranded/Coax Archwires

	ARCH FORM	Round	Rectangular	pc/Pack
TWISTADENT 3-strand <b>T3</b>	NATURAL			10
	STRAIGHT 14 inch			10
	Spool 30 foot			1
TWISTADENT 6-Coax <b>T6</b>	NATURAL			10
	STRAIGHT 14 inch			10
	Spool 30 foot			1
TWISTADENT 7-strand <b>T7</b>	Spool 30 foot			1
TWISTADENT 8-braid <b>T8</b>	NATURAL			10
	EURO			
	STRAIGHT 14 inch			10



#### Arch Forms: TWISTADENT 3-strand

Round and Spools	.0155	.0175	.0195	.0215
Square	.016 x .016			
Rectangular	.016 x .022	.017 x .025	.018 x .025	.019 x .025

#### Arch Forms: TWISTADENT 6-coax

Round and Spools	.0155	.0175	.0195	.0215
------------------	-------	-------	-------	-------

#### Arch Forms: TWISTADENT 7-strand

Spools	.0155	.0175
--------	-------	-------

#### Arch Forms: TWISTADENT 8-braid

Square	.016 x .016				
Rectangular	.016 x .022	.017 x .025	.018 x .025	.019 x .025	.021 x .025

**Single Packaged Wires available for:** EURO, NATURAL

Single packaged archwires on request available.

10 individually packaged wires per pack

\*\*\* Minimum order of 10 packs per item number required \*\*\*

## ARCHWIRE ORDER INFORMATION

In order to facilitate your ordering process, you can customize your wire and create your item number by combining the following digits:

**T3 or T6 or T7 or T8 + Arch Form + Wire Diameter + Jaw Code**

		OBLIGATORY ARCH FORM		OBLIGATORY Wire Diameter		OBLIGATORY Jaw Code	
T3 or T6 or T7 or T8		EURO	E	.012"	12	Upper	U
		NATURAL	N	.014"	14	Lower	L
		STRAIGHT	ST	.016"	16	Upper & Lower	
		SPOOL	SP	...	...	Spool	
				.019" x .025"	19x25		
				Please refer to the corresponding Arch-wire Availability to see which arch diameters are available for your selection!			
						T7 + SP + 175 = T7SP175	

**EXAMPLE 1** You would like to order a **7-braid Stainless Steel** archwire on a **SPOOL** with the diameter **.0175**.

**T8 + E + 18x25 + U = T8E18x25U**

**EXAMPLE 2** You would like to order a **8-braid Stainless Steel** archwire in **EURO** Form with the diameter **.018x.025** for the **Upper Arch**.

For Single Packaged Arch Wires please use the same Order Logic as above stated. Please put a "S-" in front of the item number to indicate that you would like to have your wire selection individually packaged.